

Gel Documentation and Analysis Automated imaging



GEL IMAGING AND ANALYSIS

Automated imaging for all your applications

Syngene imaging systems are recognised world-wide as high quality, high performance instruments for the capture and analysis of fluorescent gels, chemiluminescent western blots and protein samples. With over 25 years of experience our systems are respected for their reliability and accuracy. Syngene has pioneered new innovations in gel imaging and this latest series of the acclaimed G:BOX is no exception.

The **G:BOX** range features the new GeneSys control software which takes automation to a new level.

ALL SYNGENE SYSTEMS ARE NOW 'APPLICATION DRIVEN'.

With a **G:BOX** system you simply tell it your application and GeneSys does the rest. Using an extensive database the system calculates the best configuration

for any dye or substrate you are using and automatically sets camera, lens, lighting and filters as necessary.

> This leaves you to view the final high quality image at the touch of an on-screen button.

To extend the range of applications the G:BOX uses a unique 'plug and play' lighting system featuring LED modules*. These 'clip-in' modules are available with Red, Green, Blue and IR outputs allowing you to work with an extensive array of colourimetric dyes.

Additionally, modules for working with multiplex samples are also available.*

As usual with any Syngene system, analysis software comes as standard. The GeneTools package for 1D analysis, BP/MW, quantification, spot blots and colony counting is the first choice for many researchers and the latest version is supplied with all **G:BOX** systems.

There are 6 models to choose from across the G:BOX range.



G:BOX F3





Ultra quiet cameras ranging from 1.4m up to 8.3m pixels. Some cameras have varying degrees of cooling.



Lenses

All lenses are computer controlled and are motor driven. When the system has the motor driven stage (XL $^{1.4}$, XT 4 , IR 6 and XX 8 models) the lens automatically tracks any upward or downward movement of the sample using auto-focus. The XL $^{1.4}$, XT 4 , IR 6 and XX 8 models have exceptional high performance essential for some chemiluminescence, IR or low light level fluorescence applications.



Filters

A 7 position computer controlled filter wheel is standard. There is an extensive range of emission filters for an array of applications. Check the Syngene Product and Application selector on www.syngene.com for more information.



Overhead (Epi) illumination

White light

Every **G:BOX** has an integral white LED lighting system.

Ultra violet (UV) option

A module can be fitted with either 365, 302 or 254nm UV tubes.

LED lighting option*

To extend the range of applications, the **G:BOX** uses a unique 'plug and play' lighting system. The LED lighting modules are single wavelength and simply clip into a lighting gantry.

Up to 4 modules can be plugged into each side of the darkroom, either of the same or mixed wavelengths.

There are 9 different LED modules to choose from. Standard modules have cut off filters and generally produce a wider output wavelength whereas the 'M' series is fitted with special 'band pass' filters which allows the imaging of multiplexed blots with minimal cross-talk between dyes.



Visible light (option)

For extending the transmitted light applications such as Coomassie gels, a conversion screen is available.

Visible light converter screen for all visible light applications



Adjustable stage

On the $XL^{1.4}$, XT^4 , IR^6 and XX^8 models there is a motor driven stage included as standard. This enables samples to be moved closer to the camera. This is especially useful when working with smaller sample sizes. The stage is also linked to the GeneSys controller so that the lens automatically 'tracks' the sample as it is moved up or down providing auto-focus.



Darkroom

G:BOX has a fully light tight, robust darkroom suitable for advanced chemiluminescence, fluorescence and visible light applications. A hood is supplied with all models (optional on the **G:BOX** F³). An electronic, auto-door lock with security function prevents interruption to long exposures.



Transilluminator (option)

A 302nm transilluminator slides out for easy access. It has a variable intensity setting and has a safety cut-off. Other wavelengths are available.



UltraBright - LED transilluminator (option)

LED transilluminators operating at 470nm are suitable for use with a range of dyes including all the 'safe' dyes which are alternatives to ethidium bromide.

G:BOX APPLICATIONS

Applications you can image with a G:BOX

DNA detection with Ethidium Bromide or 'Safe' dyes

Using a 302nm UV transilluminator and Orange emission filter, images of ethidium bromide stained DNA gels can be captured in a fraction of a second. An alternative using less harmful 'Safe' dyes and blue light can produce similar images of DNA gels without loss of sensitivity.

Protein analysis

Protein gels stained with either Coomassie blue stain or silver stain can easily be captured with a **G:BOX** using a visible light converter screen with the UV transilluminator or a white light table.

Fluorescence stained gels

Syngene's LED lighting modules (Red, Green, Blue and IR) cover a wide range of fluorescence applications and dyes eg, Q-Dots, Dy lights, Alexa Fluor, CyDye and LI-COR IR dyes. Using the correct excitation light source and corresponding emission filter virtually any application can be used with a **G:BOX**. The Syngene database which controls a **G:BOX** already has thousands of different imaging protocols which it automatically selects from when using these fluorescent dyes.

Chemiluminescence with or without colourimetric markers*

Chemiluminescence with a **G:BOX** is very easy. Simply place the blot inside the system and GeneSys will assess the output. Within seconds you can make a choice between having a speedy image or one taking a little more time if high quality is required. Either way you can capture quality chemiluminescence images with accuracy and minimal effort. Colourimetric markers can also be automatically captured and merged on to the same image.

Multiplex imaging*


The use of multiplex dyes on gels or blots is increasing. With a **G:BOX** using GeneSys the user need only select the dyes being used and the system calculates the rest. Dye spectra, lighting and filters are all contained within the database. The result is that effortlessly you get perfect multiplexed images. The system uses the LED lighting modules that are configured for multiplexed dyes in conjunction with the corresponding filters.

Note: Not all **G:BOX** systems can image all of the above applications – please check the full specification before ordering.
*except **G:BOX** F³

Safe dye O
Coomassie O
Silver O

Single colour dye O
Chemi image with marker O

Multiplex imaging



POWERED BY GENESYS

Applications at the touch of a button

The biggest difference between a Syngene system and any other is in the way in which the **G:BOX** is controlled.

GENESYS, AT THE HEART OF EVERY G:BOX, IS APPLICATION DRIVEN SOFTWARE.

You will know exactly what your application is and how you have prepared your gel or blot. You simply enter this information into GeneSys and the system takes over the rest. At the heart of GeneSys is an extensive database containing data relating to a very wide range of applications, eg, fluorescence, chemiluminescence and chemifluorescence. Once the **G:BOX** has been told what sample to expect, GeneSys decides what hardware configuration is best and sets the system ready for image capture. Camera control, exposure time, sensitivity setting, lighting requirements, lens control, filter selection - all this is controlled by GeneSys.

The fully automatic functions of GeneSys allow even the most difficult or complicated gels and blots to be imaged. With chemiluminescence, the system uses its intelligent database to know the output profile of most of the major substrates. It then decides the best exposure setting. During the image capture process it continually monitors the output from the blot until the best image is produced.

GeneSys is ideal for multiplexed gels and blots. Once the application information is selected by the user, GeneSys will make separate exposures using the correct filter, lighting and camera setting for each dye. The final multiplex image will be displayed.


GeneSys has been designed for use with a touch panel PC and to take advantage of the new functions afforded by Windows 7.

GeneSys lets you control the system automatically or manually. An image browser facility allows you to view stored images. Previously saved protocol configurations can easily be accessed.

GeneSys is icon driven and has been designed for use with a touch screen although using a keyboard and mouse with a monitor is still possible.










The type of dye being used can be selected from the extensive internal database.



When using automatic capture the system will use the optimum lighting, filter and exposure time to produce the best image. Captured images are shown in the main part of the screen while previously captured images are displayed in an image pool as a series of thumbnails. The capture of multiplex images is very easy in auto-mode. The system will capture each dye separately using the optimum settings ensuring minimal crossover before displaying each image and the composite image.









The user can have full manual control of all functions.






GEL AND BLOT ANALYSIS WITH GENETOOLS

Every G:BOX comes with copies of our acclaimed GeneTools analysis software. This highly automated software can rapidly analyse a gel from loading to output of results in a matter of seconds. Simplicity is the key feature of GeneTools allowing even the most inexperienced user to obtain fast, accurate results at the click of a button and with minimal training.









- 1D gel analysis MW/BP calculation Quantification E-gels
- O Band matching with dendrograms O Spot blots O Slot blots O Colony counting
- GeneDirectory (option) for extended band matching, cluster analysis, VNTR analysis, genotyping, RFLP studies, dendrogram generation, bootstrapping



G:BOX Specifications

	0.70	0.70	0.701/	0.701/		
	G:BOX F ³	G:BOX CHEMI XR ⁵	G:BOX CHEMI XL ^{1,4}	G:BOX CHEMI XT ⁴	G:BOX CHEMI XX ⁸	G:BOX CHEMI IR6
Pixels (Million)	3.8	5.0	1.4	4.2	8.3	6.3
A/D	12/16 bit	16 bit	16 bit	16 bit	16 bit	16 bit
Greyscales	4096/65536	65536	65536	65536	65536	65536
Dynamic range	3.6/4.8	4.8	4.8	4.8	4.8	4.8
Cooling		Ultra peltier cooling for extra long exposures REGULATED	Ultra peltier cooling for extra long exposures REGULATED	Ultra 2 stage peltier cooling for extra long exposures REGULATED	Ultra peltier cooling for extra long exposures REGULATED	Ultra peltier cooling for extra long exposures REGULATED
Lens	F1.2	F1.2	F0.85	F0.95	F1.4	F1.4
Filter wheel motor driven	Yes	Yes	Yes	Yes	Yes	Yes
Illumination						
UV transilluminator (20x20cm)	V transilluminator (20x20cm) Option Op		Option	Option	Option	Option
Visible light converter Option		Option	Option	Option	Option	Option
White Epi - overhead	Yes	Yes	Yes	Yes	Yes	Yes
Epi UV	Option	Option	Option	Option	Option	Option
Epi LEDs for fluorescence	No	Option	Option	Option	Option	Option
Epi IR	No	Option	Option	Option	Option	Yes
Fold down visible light table	Option	Option	No	No	No	No
UltraBright - LED transilluminator	Option	Option	Option	Option	Option	Option
Dimensions (W,H,D cm)	57 x 85 x 43	57 x 85 x 43	57 x 100 x 52	57 × 100 × 52	57 × 100 × 52	57 × 100 × 52
Weight	37kg	37kg	45kg	45kg	45kg	45kg



G:BOX models: There are 6 G:BOX systems covering a wide range of applications

G:BOX F ³	This system features a 3.8m pixel camera and is the ideal choice for general fluorescence work such as ethidium bromide gels and proteins.
G:BOX XR ⁵	For general fluorescence and chemiluminescence applications this model with a 5.0m pixel camera is perfect. It is suitable for upgrade by adding lighting modules.
G:BOX XL ^{1,4}	The XL $^{1.4}$ has been configured for high performance chemiluminescence and uses an exceptionally sensitive 1.4m pixel camera and a large aperture (f0.85) lens. The ultra sensitive sensor has pixels of size 6.45 x 6.45um. The system can be upgraded for fluorescence applications if required.
G:BOX XT ⁴	This high performing system is suitable for fluorescence and chemiluminescence using a 4.2m pixel camera. It has a large aperture lens (f0.95) making it perfect for low light level applications where longer exposures are needed.
G:BOX XX ⁸	Our highest resolution system at 8.3m pixels for fluorescence and chemiluminescence.
G:BOX IR ⁶	Highly sensitive system suitable for IR fluorescence, fluorescence and chemiluminescence. This 6.3m pixel camera is ideal for low light level applications.

	G:BOX F³	G:BOX CHEMI XR ⁵	G:BOX CHEMI XL ^{1,4}	G:BOX CHEMI XT ⁴	G:BOX CHEMI XX8	G:BOX CHEMI IR ⁶
Chemiluminescence eg. ChemiFast, WestDura, ECL, ECL Plus	No	Yes	Yes	Yes	Yes	Yes
Fluorescence Using transmitted 302nm UV Ethidium Bromide, Gel Red, ProQ Diamond, Deep Purple	Yes	Yes	Yes	Yes	Yes	Yes
Lighting options						
Fluorescence Using Epi 365nm UV option QDots, TLC, ECL Plus	Yes	Yes	Yes	Yes	Yes	Yes
Fluorescence Using Epi Blue LED (480nm) option Alexafluor 488, Cy2	No	Yes	Yes	Yes	Yes	Yes
Fluorescence Using Epi Green LED (540nm) option Alexafluor 546, Cy3, HEX, TAMRA	No	Yes	Yes	Yes	Yes	Yes
Fluorescence Using Epi Red LED (640nm) option Alexafluor 647, Alexafluor 633, Cy5, BODIPY 650	No	Yes	Yes	Yes	Yes	Yes
Fluorescence Using Epi Near IR (740nm) option Alexafluor 700, Alexafluor 750, Cy7	No	Yes	Yes	Yes	Yes	Yes
Fluorescence Using Epi-red LED (640nm) option IR Dye 680, IR Dye 700	No	Yes	Yes	Yes	Yes	Yes
Fluorescence Using Epi IR (740nm) option IR Dye 800, AlexaFluor 790	No	Yes	Yes	Yes	Yes	Yes
Visible light imaging Using conversion screen option with transilluminator Coomassie blue, Silver stain, X-ray film	Yes	Yes	Yes	Yes	Yes	Yes
Blue light imaging Using UltraBright - LED blue transilluminator option SYBR stains, Gel Green, UltraSafe blue	Yes	Yes	Yes	Yes	Yes	Yes

Over 75,000 scientists world-wide in pharmaceutical and biotech companies, as well as academic and government institutions, have chosen Syngene as their expert imaging partner. If you'd like to find out why, please contact us or one of our dealers for more information and a demonstration of the revolutionary **G:BOX**

Syngene Europe and International Headquarters:

Beacon House Nuffield Road Cambridge CB4 1TF UK Tel: +44 (0)1223 727123 Fax: +44 (0)1223 727101 email: sales@syngene.com

Syngene USA Headquarters:

5108 Pegasus Court Suite M Frederick MD 21704 USA Tel: 800-686-4407/301-662-2863 Fax: 301-631-3977

email: ussales@syngene.com

Website: www.syngene.com



Please refer to www.syngene.com for all ordering information

