

Cellometer[®] Auto 1000

Automated Cell Viability Counter

Includes:

- Cellometer Auto 1000 Instrument with integrated computer and touch-screen
- Cellometer Software
- USB 2.0 Connection Cable
- Power Supply
- Ethernet Cable
- Phone/online applications support during set-up

**Imaging
Performance:**

Cell Size 5 - 60 microns (most mammalian cells)
Conc. Range: 10⁵ – 10⁷ cells/ml

**Instrument
Specifications:**

Weight: 9 lbs. (4.1 kg)

Dimensions: Width: 11.1" (28.2 cm)
 Depth: 8.3" (21.1 cm)
 Height: 9.0" (22.9 cm)

Electrical Input

Power Adapter: 100-240 VAC, 50-60 Hz, 1.5A
Input to Instrument: 12 VDC, 7.5A

Display:

10.4" touch screen

For more information, email support@nexcelom.com or call: 978-327-5340.



Nexcelom Bioscience LLC. | 360 Merrimack Street, Building 9 | Lawrence, MA 01843
Telephone: 978.327.5340 | Fax: 978.327.5341 | Email: info@nexcelom.com | www.nexcelom.com