

Cellometer® Auto 1000

Cell Counter

Accurate Cellometer Cell Counting ... Now Even Faster and Easier

- ➔ Improve counting of live/dead cells, individual cells within clusters, and irregular-shaped cells: Proprietary Cellometer software
- ➔ Spend less time counting and more time on research: Concentration and viability in <10 seconds!
- ➔ Easily count cells in the hood: Stand-alone system with a large, user-friendly touch screen
- ➔ View, print, and save counted cell images, data reports, and cell size histograms: Connect to a network or use a USB drive



**U.S. Introductory
Price of \$4,998**
valid through
April 30, 2013.

Simple, Automated Cell Counting


Improve Speed, Accuracy, and Consistency



Load 20 µl of Sample




Insert Counting Chamber



Select Assay & Click Count



View Image & Results



For more information on the new Cellometer Auto 1000, contact Nexcelom at 978-327-5340 or info@nexcelom.com!